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Referências

ACC/ AHA/ ESC 2006 Guidelines for the Management of Patients with Atrial Fibrillation. *Circulation* 2006;114:700-752.

Lown, B. et al. "Cardioversion" of atrial fibrillation. *New England Journal of Medicine* 1963;269:325-331.

Field, J.M., Hazinski, M.F., & Gilmore, D. [Eds]. *2006 Handbook of Emergency Cardiovascular Care for Healthcare Providers*. American Heart Association.

Neal, S., Ngarmukos, T., Lessard, D. & Rosenthal, L. Comparison of the efficacy and safety of two biphasic defibrillator waveforms for the conversion of atrial fibrillation to sinus rhythm. *American Journal of Cardiology* 2003;92:810-814.

Mittal, S Ayati, S et al. Transthoracic Cardioversion of Atrial Fibrillation. Comparison of Rectilinear Biphasic Versus Damped Sine Wave Monophasic Shocks. *Circulation*. 2000;101:1282-1287

Niebauer M et al. Comparison of the rectilinear biphasic waveform with the damped sine monophasic waveform for external cardioversion of atrial fibrillation and flutter. *Am J Card*. 2004: Vol. 93. 1495 -99

Sado, D.M. et al. Comparison of the effects of removal of chest hair with not doing so before external defibrillation on transthoracic impedance. *American Journal of Cardiology* 2004;93:98-100.

Sirna, S.I et al. Factors affecting transthoracic impedance during electrical cardioversion. *American Journal of Cardiology* 1988;62:1048-1052.

Kirchhof, P. et al. Anterior-posterior versus anterior-lateral electrode positions for external cardioversion of atrial fibrillation: A randomized trial. *Lancet* 2002;360:1275-1279.

Sirna, S.I et al. Factors affecting transthoracic impedance during electrical cardioversion. *American Journal of Cardiology* 1988;62:1048-1052.

Mehdirad AA et al Improved Clinical Efficacy of External Cardioversion by Fluoroscopic Electrode Positioning and Comparison to Internal Cardioversion in Patients with Atrial Fibrillation *PACE* 1999 Volume 22 Issue 1, 237